


PILLAR HMA I-130
45016-XXXXX-XXXXXXX



Stable pillar for rigid mounting of HMA I-130 suspension system on floor or ceiling


- Default colours
 - Components: RAL 7035 (Light grey)
 - Special colours on request
- Internal cable ducting

Product code 45016-xxxxx-xxxxxxx

Safety Class IP 54 DIN EN 60529

FLOOR PILLAR HMA I-130

45016-XXXXX-XXXXXXX



Stable pillars for rigid mounting of HMA I-130 suspension systems on floor or ceiling

- Default colour
 - RAL 7035 (Light grey)
 - Special colours on request
- Height up to 2500 mm in 100 mm increments
- Internal cable ducting
- Top flange plate for HMA I-130 top mounted joint

Product code	Description
45016-05000-1703511	Floor pillar HMA I-130, 500 , Colour: RAL 7035
45016-06000-1703511	Floor pillar HMA I-130, 600 , Colour: RAL 7035
45016-07000-1703511	Floor pillar HMA I-130, 700 , Colour: RAL 7035
45016-08000-1703511	Floor pillar HMA I-130, 800 , Colour: RAL 7035
45016-09000-1703511	Floor pillar HMA I-130, 900 , Colour: RAL 7035
45016-10000-1703511	Floor pillar HMA I-130, 1000 , Colour: RAL 7035
45016-11000-1703511	Floor pillar HMA I-130, 1100 , Colour: RAL 7035
45016-12000-1703511	Floor pillar HMA I-130, 1200 , Colour: RAL 7035
45016-13000-1703511	Floor pillar HMA I-130, 1300 , Colour: RAL 7035
45016-14000-1703511	Floor pillar HMA I-130, 1400 , Colour: RAL 7035
45016-15000-1703511	Floor pillar HMA I-130, 1500 , Colour: RAL 7035
45016-16000-1703511	Floor pillar HMA I-130, 1600 , Colour: RAL 7035
45016-17000-1703511	Floor pillar HMA I-130, 1700 , Colour: RAL 7035
45016-18000-1703511	Floor pillar HMA I-130, 1800 , Colour: RAL 7035
45016-19000-1703511	Floor pillar HMA I-130, 1900 , Colour: RAL 7035
45016-20000-1703511	Floor pillar HMA I-130, 2000 , Colour: RAL 7035
45016-21000-1703511	Floor pillar HMA I-130, 2100 , Colour: RAL 7035
45016-22000-1703511	Floor pillar HMA I-130, 2200 , Colour: RAL 7035
45016-23000-1703511	Floor pillar HMA I-130, 2300 , Colour: RAL 7035
45016-24000-1703511	Floor pillar HMA I-130, 2400 , Colour: RAL 7035
45016-25000-1703511	Floor pillar HMA I-130, 2500 , Colour: RAL 7035