



HMA LIFT 20



Continuous height adjustment suspension arm for medical equipment like microscopes, control panel, etc.

- Default colours
 - Cover: RAL 9016 (Traffic white)
 - Special colours on request
- Adjustable stroke limit
- Stroke locking with operator release (option)
- Ease of vertical movement separately adjustable (option)
- Restraint-guided load suspension at head joint
- Fine adjustment of the load capacity range (TLB) using adjusting screw
- external cable ducting system
- Limiting of the swivelling range (option)
- Adjustable rotation limit of the horizontal swivel angle (option)
- Head bearing mounting upwards and downwards possible
- Low maintenance

Cable duct	In external cable duct (Internal cable duct 8,2 x 14mm option)
Vertical stroke	0.00
Loading weight range	1.5 - 34 kg

CABLE COVERING LIFT (360 MM)

4-6-25-801



Cable covering for lift units in four different lengths

- Default colour
 - RAL 7035 (Light grey)
 - Special colours on request

HMA I-020 KEYBOARD TRAY

62581-20070-XXXXXXX



For mounting on an extension arm or a lift unit

- Default colour
 - RAL 7035 (Light grey)
 - Special colours on request

Product code	Description
62581-20070-1703510	HMA I-020 Keyboard Tray, 4kg, Colour: RAL 7035

HMA I-020 KEYBOARD TRAY WITH MOUSEPAD

62581-20080-XXXXXXX



For mounting on an extension arm or a lift unit

- Default colour
 - RAL 7035 (Light grey)
 - Special colours on request

Product code	Description
62581-20080-1703510	HMA I-020 Keyboard Tray with Mousepad, 4.7kg, Colour: RAL 7035

HMA I-020 SWIVEL- / TILT JOINT L

62581-10030-XXXXXXX



For mounting on an extension arm or unit

- Default colour
 - RAL 7035 (Light grey)
 - Special colours on request
- Swing angle: 360°
- Swivel range: up to 180°
- Connection according VESA®75, optional VESA®100 possible

Product code	Description
62581-10030-1703510	HMA I-020 Swivel- / Tilt Joint L, 0.9kg, Colour: RAL 7035

HMA I-020 SWIVEL- / TILT JOINT M

62581-10020-XXXXXXX



For mounting on an extension arm or unit

- Default colour
 - RAL 7035 (Light grey)
 - Special colours on request
- Swing angle: 360°
- Swivel range: up to 180°
- Connection according VESA® 75 or VESA® 100

Product code	Description
62581-10020-1703510	HMA I-020 Swivel- / Tilt Joint M, 0.85kg, Colour: RAL 7035

HMA I-020 SWIVEL- / TILT JOINT M-R

62581-10040-XXXXXXX



For mounting on an extension arm or unit

- Default colour
 - RAL 7035 (Light grey)
 - Special colours on request
- Angle of rotation: 360°
- Rotation: 360°
- Swivel range: up to 180°
- Connection according VESA® 75 or VESA® 100

Product code	Description
62581-10040-1703510	HMA I-020 Swivel- / Tilt Joint M-R, 1kg, Colour: RAL 7035

SWIVEL- / TILT JOINT S (81 MM)

4-6-25-601



For mounting on an extension arm or unit

- Default colour
 - RAL 7035 (Light grey)
 - Special colours on request
- Swing angle: 360°
- Swivel range: up to 180°
- Connection according VESA® 75 or VESA® 100

SWIVEL- / TILT JOINT S-R (103 MM)

4-6-25-605



For mounting on an extension arm or unit

- Default colour
 - RAL 7035 (Light grey)
 - Special colours on request
- Angle of rotation: 360°
- Rotation: 360°
- Swivel range: up to 180°
- Connection according VESA® 75 or VESA® 100

HMA I-020 ANGLE ADAPTOR

62581-40020-XXXXXXX



For rear mounting of an enclosure on a lift unit

- Default colour
 - RAL 7035 (Light grey)
 - Special colours on request

Product code	Description
62581-40020-1703510	HMA I-020 Angle Adaptor, 1.35kg, Colour: RAL 7035

HMA I-020 EXTENSION ARM

62571-110XX-XXXXXXX



Extension arm in 9 different lengths

- Default colour
 - RAL 7035 (Light grey)
 - Special colours on request
- Overhead-mounting possible

Product code	Description
62571-11012-1703510	HMA I-020 Extension Arm, 120 mm, 0.45kg, Colour: RAL 7035
62571-11015-1703510	HMA I-020 Extension Arm, 150 mm, 0.5kg, Colour: RAL 7035
62571-11017-1703510	HMA I-020 Extension Arm, 170 mm, 0.55kg, Colour: RAL 7035

Product code	Description
62571-11019-1703510	HMA I-020 Extension Arm, 190 mm, 0.6kg, Colour: RAL 7035
62571-11025-1703510	HMA I-020 Extension Arm, 250 mm, 0.7kg, Colour: RAL 7035

HMA I-020 ADAPTOR

62581-30010-XXXXXXX



For mounting an enclosure

- Default colour
 - RAL 7035 (Light grey)
 - Special colours on request

Product code	Description
62581-30010-1703510	HMA I-020 Adaptor, 0.55kg, Colour: RAL 7035