



HMA M40



- Default colours
 - Pipes: RAL 9016 (Traffic white)
 - Components: RAL 9016 (Traffic white)
 - Special colours on request

Safety Class

IP 54 DIN EN 60529

CEILING MOUNTING



For mounting lift systems on the ceiling

- Default colour
 - RAL 9016 (Traffic white)
 - Special colours on request

U-EXTENSION ARM



System extension arm for wall mounting

- Default colour
 - RAL 9016 (Traffic white)
 - Special colours on request

WALL MOUNTING FOR STANDING MOUNTING

500-02-10100-9016



For hanging or standing wall mounting

- Default colour
 - RAL 9016 (Traffic white)
 - Special colours on request

WALL MOUNTING FOR HANGING MOUNTING

500-02-10101-9016



For hanging or standing wall mounting

- Default colour
 - RAL 9016 (Traffic white)
 - Special colours on request

Z-EXTENSION ARM



System extension arm for ceiling mounting

- Default colour
 - RAL 9016 (Traffic white)
 - Special colours on request

INCLINATION ADAPTOR (144 MM)

500-05-00101-9016



For movable connection of attachment parts with the lift

- Default colour
 - RAL 9016 (Traffic white)
 - Special colours on request
- Vertical swivel range $\pm 90^\circ$
- Dimensions of the flange panel: 60 x 90 mm

INCLINATION ADAPTOR VESA®

500-05-00110-9016



For adjustable connection of monitors with VESA® connection

- Default colour
 - RAL 9016 (Traffic white)
 - Special colours on request
- Vertical swivel range: 25°
- Removable lid for easy installation

TUBE ADAPTOR (300 MM)

500-05-00200-0000



For connection of a circular tube Ø 25 mm

- Default colour
 - Stainless steel
 - Special colours on request